

Anti-Argonaute-2 (Eif2C2) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukv8

Item Number: CM0030687



Introduction

High-quality Rabbit Polyclonal Antibody against human Argonaute-2 (AGO2/EIF2C2). Validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for RNAi and gene silencing research.

[Learn More](#)

Parameter	Value
Product Name	Argonaute-2 Rabbit pAb
Remarks/Alias	PPD; Q10; CASC7; EIF2C2; LESKRES; LINC00980; Argonaute-2
Species	Human
Gene ID (Human)	27161
Gene ID	27161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-115 of human Argonaute-2 (NP_036286.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKV8
Protein Weight	97kDa
Shipping	Ice bag