

Anti-Argonaute-2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ukv8

Item Number: CM0002445



Introduction

High-quality rabbit monoclonal antibody targeting human Argonaute-2 (Ago2). Validated for Western blot, immunohistochemistry, and ELISA with cross-reactivity to mouse and rat. Ideal for RNAi and gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	Argonaute-2 Rabbit mAb
Remarks/Alias	PPD; Q10; CASC7; EIF2C2; LESKRES; LINC00980; Ago 2; Argonaute-2
Species	Human
Gene ID (Human)	27161
Gene ID	27161
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKV8
Protein Weight	97kDa
Shipping	Ice bag