


Anti-Arid3A Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q99856

Item Number: CM0026935



Introduction

ARID3A rabbit polyclonal antibody for WB, IHC-P, IP, ELISA. Reacts with human, mouse, rat. Validated for B-cell differentiation and cell cycle research. Protein weight 63kDa.

[Learn More](#)

Parameter	Value
Product Name	ARID3A Rabbit pAb
Remarks/Alias	DRIL1; DRIL3; BRIGHT; E2FBP1; ARID3A
Species	Human
Gene ID	1820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-593 of human ARID3A (NP_005215.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99856
Protein Weight	63kDa
Shipping	Ice bag