


Anti-Arp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61158

Item Number: CM0013957



Introduction

Rabbit polyclonal antibody targeting human Arp3 (ACTR3), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse. For actin polymerization, cell motility, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Arp3 Rabbit pAb
Remarks/Alias	ARP3; Arp3
Species	Human
Gene ID	10096 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-418 of human Arp3 (NP_005712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61158
Protein Weight	47kDa
Shipping	Ice bag