

Anti-Arsk Polyclonal Antibody For Wb, Elisa - Q9D2L1

Item Number: CM0022533

Introduction

Arsk Rabbit pAb validated for WB and ELISA applications. Targets human arylsulfatase K with recombinant immunogen. Ideal for lysosomal enzyme and glycosaminoglycan degradation studies.

[Learn More](#)

Parameter	Value
Product Name	Arsk Rabbit pAb
Remarks/Alias	ASK; 2810429K17Rik; 4833414G15Rik; Arsk
Species	Human
GeneID (human)	77041
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-553 of human Arsk(NP_084123.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9D2L1
Protein Weight	63kDa
Shipping	Ice bag