

Anti-Asb2 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Q27

Item Number: CM0022811

Introduction

High-quality rabbit polyclonal antibody targeting human ASB2, validated for WB and ELISA with reactivity in human and rat. Ideal for studying ubiquitin-mediated degradation in cardiac and hematopoietic research. SWISS: Q96Q27.

[Learn More](#)

Parameter	Value
Product Name	ASB2 Rabbit pAb
Remarks/Alias	ASB-2; ASB2
Species	Human
Gene ID (Human)	51676
Gene ID	51676
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 448-587 of human ASB2 (NP_057234.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96Q27
Protein Weight	70kDa
Shipping	Ice bag