


Anti-Asgr1 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P07306

Item Number: CM0014510



Introduction

High-quality rabbit polyclonal antibody targeting human ASGR1, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying endocytosis and SARS-CoV-2 entry.

[Learn More](#)

Parameter	Value
Product Name	ASGR1 Rabbit pAb
Remarks/Alias	HL-1; ASGPR; ASGPR1; CLEC4H1; ASGR1
Species	Human
GeneID (Human)	432
GeneID	432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-291 of human ASGR1 (NP_001662.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07306
Protein Weight	33kDa
Shipping	Ice bag