

Anti-Atf4 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P18848

Item Number: CM0019521



Introduction

Rabbit polyclonal antibody targeting human ATF4 (CREB2). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. ATF4 is a key transcription factor in cellular stress response and metabolism. Ideal for ER stress, oxidative stress, and bone development research.

[Learn More](#)

Parameter	Value
Product Name	ATF4 Rabbit pAb
Remarks/Alias	CREB2; TXREB; CREB-2; TAXREB67; ATF4
Species	Human
GeneID (Human)	468
GeneID	468
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 260-351 of human ATF4 (NP_001666.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18848
Protein Weight	39kDa
Shipping	Ice bag