

Anti-Atf5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, And Elisa - Q9Y2D1

Item Number: CM0021088

Introduction

Rabbit polyclonal antibody against human ATF5, validated for WB, IHC-P, IF/ICC, and ELISA. Detects a band at ~31 kDa in Western blot. Reacts with human, mouse, and rat. Suitable for transcription regulation and neurobiology research.

[Learn More](#)

Parameter	Specification
Product Name	ATF5 Rabbit pAb
Remarks/Alias	ATFX; HMFN0395; ATF5
Species	Human
Gene ID (Human)	22809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-282 of human ATF5 (NP_036200.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2D1
Protein Weight	31kDa
Shipping	Ice bag