

Anti-Atg16L2 Polyclonal Antibody For Wb And Elisa - Q8Naa4

Item Number: CM0014942

Introduction

ATG16L2 Rabbit Polyclonal Antibody validated for Western blot and ELISA, reacting with human, mouse, and rat samples. Suitable for autophagy and epithelial homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	ATG16L2 Rabbit pAb
Remarks/Alias	WDR80; ATG16B; ATG16L2
Species	Human
GeneID (Human)	89849
GeneID	89849
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human ATG16L2 (NP_203746.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NAA4
Protein Weight	69kDa
Shipping	Ice bag