

Anti-Atoh7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N100

Item Number: CM0030905

Introduction

High-quality rabbit polyclonal antibody targeting human ATOH7 (bHLHa13), suitable for WB, IF/ICC, and ELISA applications. Ideal for investigating retinal ganglion cell development and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	ATOH7 Rabbit pAb
Remarks/Alias	Math5; NCRNA; RNANC; PHPVAR; bHLHa13; ATOH7
Species	Human
GeneID	220202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-152 of human ATOH7 (NP_660161.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8N100
Protein Weight	17kDa
Shipping	Ice bag