

Anti-Atoh8 Polyclonal Antibody For Wb And Elisa - Q96Sq7

Item Number: CM0024321

Introduction

Rabbit polyclonal antibody targeting human ATOH8, a bHLH transcription factor. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studying cell differentiation, endothelial function, and developmental processes.

[Learn More](#)

Parameter	Specification
Product Name	ATOH8 Rabbit pAb
Remarks/Alias	HATH6; bHLHa21; ATOH8
Species	Human
Gene ID (Human)	84913
Gene ID	84913
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human ATOH8 (NP_116216.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96SQ7
Protein Weight	35kDa
Shipping	Ice bag