


Anti-Atp13A2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nq11

Item Number: CM0015727



Introduction

ATP13A2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat ATP13A2 (~129 kDa). Suitable for Parkinson's disease and lysosomal function research.

[Learn More](#)

Parameter	Value
Product Name	ATP13A2 Rabbit pAb
Remarks/Alias	CLN12; KRPPD; PARK9; SPG78; HSA9947; ATP13A2
Species	Human
GeneID (human)	23400
GeneID	23400
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 70-210 of human ATP13A2 (NP_071372.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ11
Protein Weight	129kDa
Shipping	Ice bag