

Anti-Atp1A2 Monoclonal Antibody For Wb, Elisa - P50993

Item Number: CM0003051

Introduction

Rabbit monoclonal antibody targeting human ATP1A2 (sodium/potassium-transporting ATPase subunit alpha-2). Validated for WB and ELISA with cross-reactivity to mouse and rat. Research use only.

[Learn More](#)

Parameter	Value
Product Name	ATP1A2 Rabbit mAb
Remarks/Alias	FHM2; MHP2; DEE98; FARIMPD; ATP1A2
Species	Human
GeneID (human)	477
GeneID	477
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P50993
Protein Weight	112kDa
Shipping	Ice bag