

Anti-Atp2B4 Polyclonal Antibody For Wb, Elisa - P23634

Item Number: CM0015087



Introduction

Anti-ATP2B4 rabbit polyclonal antibody validated for WB and ELISA. Targets plasma membrane calcium ATPase 4 (PMCA4, P23634) for cell signaling and sperm motility research with mouse cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	ATP2B4 Rabbit pAb
Remarks/Alias	MXRA1; PMCA4; ATP2B2; PMCA4b; PMCA4x; ATP2B4
Species	Human
GeneID (Human)	493
GeneID	493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1141-1205 of human ATP2B4 (NP_001675.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P23634
Protein Weight	138kDa
Shipping	Ice bag