

Anti-Atp4B Rabbit Polyclonal Antibody (Gastric H⁺/K⁺ Atpase B Subunit) For Wb And Elisa - P51164

Item Number: CM0013852

Introduction

High-quality rabbit polyclonal antibody against ATP4B, the gastric H⁺/K⁺ ATPase beta subunit (33 kDa). Validated for WB and ELISA across human, mouse, and rat. Ideal for studying gastric acid secretion and parietal cell function.

[Learn More](#)

Parameter	Value
Product Name	ATP4B Rabbit pAb
Remarks/Alias	ATP6B; ATP4B
Species	Human
GeneID (Human)	496
GeneID	496
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-291 of human ATP4B (NP_000696.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51164
Protein Weight	33kDa
Shipping	Ice bag