

Anti-Atp5G1 Monoclonal Antibody For Wb, If-P, Elisa - P05496

Item Number: CM0006674



Introduction

High-quality rabbit monoclonal antibody targeting human ATP5G1, a mitochondrial ATP synthase subunit C1. Suitable for WB, IF-P, and ELISA applications with cross-reactivity in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	ATP5G1 Rabbit mAb
Remarks/Alias	ATP5A; ATP5G; ATP5G1
Species	Human
GeneID (Human)	516
GeneID	516
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATP5G1 (P05496).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05496
Protein Weight	14kDa
Shipping	Ice bag