

Anti-Atp5G1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05496

Item Number: CM0002782



Introduction

Rabbit monoclonal antibody targeting human ATP5G1, also known as ATP synthase subunit c, recommended for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for mitochondrial bioenergetics studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5G1 Rabbit mAb
Remarks/Alias	ATP5A; ATP5G; ATP5G1
Species	Human
GeneID	516
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 57-136 of human ATP5G1 (P05496)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05496
Protein Weight	14kDa
Shipping	Ice bag