


Anti-Atp5J Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18859

Item Number: CM0026588



Introduction

High-quality ATP5J rabbit polyclonal antibody for WB, IHC-P, IF/ICC, and ELISA applications. Recognizes ATP synthase peripheral stalk subunit F6 in human, mouse, and rat samples, supporting mitochondrial and energy metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ATP5J Rabbit pAb
Remarks/Alias	F6; CF6; ATP5; ATPM; ATP5A; ATP5J
Species	Human
Gene ID (Human)	522
Gene ID	522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human ATP5J (NP_001676.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18859
Protein Weight	13kDa
Shipping	Ice bag