

Anti-Atp6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P00848

Item Number: CM0030522

Introduction

Rabbit polyclonal antibody targeting ATP6 (ATP synthase F(0) complex subunit a, also known as MT-ATP6) for Western blot (WB), immunofluorescence (IF/ICC), immunohistochemistry (IHC-P), and ELISA applications. Detects human and mouse ATP6.

[Learn More](#)

Parameter	Value
Product Name	ATP6 Rabbit pAb
Remarks/Alias	MT-ATP6, MTATP6, ATPASE6
Species	Mouse
Gene ID (Human)	4508
Gene ID (Mouse)	17705
Immunogen	Synthetic peptide corresponding to amino acids 127-226 of mouse ATP6 (NP_904333.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00848
Protein Weight	25kDa
Shipping	Ice bag