


Anti-Atpb Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06576

Item Number: CM0002982



Introduction

High-affinity rabbit monoclonal antibody targeting human ATP synthase beta subunit (ATPB) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial energy metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ATPB Rabbit mAb
Remarks/Alias	ATP5B; ATPMB; ATP5B; HUMOP2; HEL-S-271; ATPB
Species	Human
GeneID (Human)	506
GeneID	506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-529 of human ATPB (NP_001677.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06576
Protein Weight	57kDa
Shipping	Ice bag