


Anti-Atrip Polyclonal Antibody For Wb, Elisa - Q8Wxe1

Item Number: CM0027773



Introduction

High-quality rabbit polyclonal antibody against human ATRIP (ATR-interacting protein) for Western blot and ELISA, with cross-reactivity to mouse. Essential for DNA damage checkpoint signaling research.

[Learn More](#)

Parameter	Value
Product Name	ATRIP Rabbit pAb
Remarks/Alias	ATRIP
Species	Human
GeneID (Human)	84126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ATRIP (NP_569055.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WXE1
Protein Weight	86kDa
Shipping	Ice bag