


Anti-Aurka Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14965

Item Number: CM0022433



Introduction

Rabbit polyclonal antibody against human AURKA (Aurora kinase A), validated for WB, IF/ICC, ELISA, cross-reacts with Human, Mouse. Key mitotic kinase involved in cell cycle regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	AURKA Rabbit pAb
Remarks/Alias	AIK; ARK1; AURA; BTAK; STK6; STK7; STK15; PPP1R47; AURKA
Species	Human
Gene ID (Human)	6790
Gene ID	6790
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 301-403 of human AURKA (NP_003591.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14965
Protein Weight	46kDa
Shipping	Ice bag