

Anti-Aurora B Rabbit Polyclonal Antibody For Wb, Elisa - Q96Gd4

Item Number: CM0016246



Introduction

CamelBio offers rabbit polyclonal antibody against Aurora B (AURKB), validated for WB and ELISA with human, mouse, rat reactivity. Essential for studying mitotic regulation and chromosomal passenger complex function.

[Learn More](#)

Parameter	Value
Product Name	Aurora B Rabbit pAb
Remarks/Alias	AIK2; AIM1; ARK2; AurB; IPL1; STK5; AIM-1; ARK-2; STK-1; STK12; PPP1R48; aurkb-sv1; aurkb-sv2; Aurora B
Species	Human
Gene ID (Human)	9212
Gene ID	9212
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aurora B (NP_004208.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GD4
Protein Weight	39kDa
Shipping	Ice bag