

Anti-Baf57/Smarce1 Rabbit Mab For Wb, If/Icc, Ip, Elisa - Q969G3

Item Number: CM0012351



Introduction

Rabbit monoclonal antibody targeting human BAF57/SMARCE1, validated for Western blot, immunofluorescence/ICC, immunoprecipitation, and ELISA. Cross-reacts with mouse and rat. Detects a chromatin remodeling subunit crucial in neurogenesis and gene regulation. Ideal for epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	BAF57/SMARCE1 Rabbit mAb
Remarks/Alias	CSS5; BAF57; BAF57/SMARCE1
Species	Human
Gene ID (Human)	6605
Gene ID	6605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-100 of human BAF57/SMARCE1 (NP_003070.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969G3
Protein Weight	47kDa
Shipping	Ice bag