

Anti-Bag2 Polyclonal Antibody For Wb & Elisa - 095816

Item Number: CM0030400

Introduction

CamelBio's Anti-BAG2 Rabbit Polyclonal Antibody targets human BAG2, a nucleotide-exchange factor for HSP70/HSC70 chaperones. Validated for Western blot and ELISA; supports proteostasis, protein quality control, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	BAG2 Rabbit pAb
Remarks/Alias	BAG-2; dj41711.2; BAG2
Species	Human
GeneID (Human)	9532
GeneID	9532
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human BAG2 (NP_004273.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	095816
Protein Weight	24kDa
Shipping	Ice bag