


Anti-Bcap29 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhq4

Item Number: CM0030881



Introduction

BCAP29 Rabbit pAb detects human, mouse, and rat BCAP29 (B-cell receptor-associated protein 29) in WB, IHC-P, and ELISA. Suitable for studying ER-Golgi transport and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	BCAP29 Rabbit pAb
Remarks/Alias	B29; BAP29; BCAP29
Species	Human
GeneID	55973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-241 of human BCAP29 (NP_061332.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHQ4
Protein Weight	28kDa
Shipping	Ice bag