

Anti-Bcat2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O15382

Item Number: CM0006857



Introduction

BCAT2 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Detects human BCAT2 and cross-reacts with mouse and rat. Ideal for branched-chain amino acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	BCAT2 Rabbit mAb
Remarks/Alias	BCAT2; BCAM; BCATM; BCT2; PP18; branched chain amino acid transaminase 2
Species	Human
Gene ID (human)	587
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-392 of human BCAT2. (NP_001181.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15382
Protein Weight	33kDa/44kDa
Shipping	Ice bag