


Anti-Bckdhb (E1B) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P21953

Item Number: CM0006536



Introduction

Rabbit monoclonal antibody against human BCKDHB (E1B, BCKDH E1-beta). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial metabolism and Maple syrup urine disease research.

[Learn More](#)

Parameter	Value
Product Name	BCKDHB Rabbit mAb
Remarks/Alias	E1B; BCKDE1B; BCKDH E1-beta; BCKDHB
Species	Human
Gene ID (Human)	594
Gene ID	594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-392 of human BCKDHB (NP_898871.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21953
Protein Weight	43 kDa
Shipping	Ice bag