

Anti-Bcl-2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P49950

Item Number: CM0010753

Introduction

Ideal for apoptosis and cancer research, this anti-Bcl-2 rabbit monoclonal antibody is validated for WB, IHC-P, and ELISA, with confirmed cross-reactivity to mouse and rat. - P49950

[Learn More](#)

Parameter	Specification
Product Name	Bcl-2 Rabbit mAb
Remarks/Alias	Bcl-2
Species	Rat
GeneID (Human)	596
GeneID	24224
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of rat Bcl2 (NP_058689.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49950
Protein Weight	26 kDa
Shipping	Ice bag