

Anti-Beclin 1 Polyclonal Antibody For Wb, Elisa - Q14457

Item Number: CM0021615



Introduction

Beclin 1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA. Detects a ~52kDa band in human, mouse, and rat samples. Suitable for autophagy and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Beclin 1 Rabbit pAb
Remarks/Alias	ATG6; VPS30; beclin1; Beclin 1
Species	Human
Gene ID (Human)	8678
Gene ID	8678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-280 of human Beclin 1 (NP_003757.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14457
Protein Weight	52kDa
Shipping	Ice bag