

Anti-Bin2 Rabbit Polyclonal Antibody For Wb And Elisa - Q39011

Item Number: CM0019531

Introduction

Affinity-purified rabbit polyclonal antibody targeting Arabidopsis thaliana BIN2 (shaggy-related protein kinase eta). Suitable for Western blot and ELISA. Negative regulator of brassinosteroid signaling. UniProt Q39011.

[Learn More](#)

Parameter	Value
Product Name	BIN2 Rabbit pAb
Remarks/Alias	ATSK21; BRASSINOSTEROID-INSENSITIVE 2; DWARF 12; DWF12; F28A21.120; F28A21_120; SHAGGY-LIKE KINASE 21; SK21; UCU1; ULTRACURVATA 1; BIN2
Species	Arabidopsis thaliana
Gene ID	827605
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 281-380 of Arabidopsis thaliana BIN2 (NP_193606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q39011
Protein Weight	43kDa
Shipping	Ice bag