


Anti-Bmp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12643

Item Number: CM0013079



Introduction

High-quality BMP2 rabbit polyclonal antibody validated for Western blot, IHC, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for osteogenesis and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	BMP2 Rabbit pAb
Remarks/Alias	BDA2; BMP2A; SSFSC; SSFSC1; BMP2
Species	Human
GeneID (Human)	650
GeneID	650
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 283-396 of human BMP2 (NP_001191.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12643
Protein Weight	45kDa
Shipping	Ice bag