


Anti-Bmpr1B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00238

Item Number: CM0023181



Introduction

Rabbit polyclonal BMPR1B antibody validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for bone morphogenetic protein signaling and chondrocyte differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	BMPR1B Rabbit pAb
Remarks/Alias	ALK6, AMD3, AMDD, BDA2, ALK-6, BDA1D, CDw293, BMPR1B
Species	Human
Gene ID	658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human BMPR1B (NP_001194.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00238
Protein Weight	57kDa
Shipping	Ice bag