


## Anti-Bms1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14692

Item Number: CM0016441



### Introduction

Rabbit polyclonal antibody targeting human BMS1 (ribosome biogenesis protein BMS1 homolog). Suitable for Western blot and ELISA; cross-reacts with mouse and rat. Detects 146 kDa protein; ideal for nucleolar and ribosome assembly studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | BMS1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | ACC; BMS1L; BMS1                                                                                                  |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 9790                                                                                                              |
| GeneID           | 9790                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human BMS1 (NP_055568.3) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q14692                                                                                                            |
| Protein Weight   | 146kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |