


Anti-Brd9 Rabbit Monoclonal Antibody For Wb, Elisa - Q9H8M2

Item Number: CM0006389



Introduction

Recombinant rabbit monoclonal antibody targeting human BRD9. Suitable for Western blot and ELISA applications; cross-reacts with mouse and rat. Target is a chromatin reader involved in transcription regulation and homologous recombination.

[Learn More](#)

Parameter	Value
Product Name	BRD9 Rabbit mAb
Remarks/Alias	PRO9856; SMARCI2; LAVS3040; BRD9
Species	Human
Gene ID (human)	65980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 274-597 of human BRD9 (NP_076413.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8M2
Protein Weight	67kDa
Shipping	Ice bag