

Anti-Btaf1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa, Chip - O14981

Item Number: CM0028153

Introduction

BTAf1 Rabbit Polyclonal Antibody (O14981) validated for WB, IP, ELISA, ChIP. Reacts with Human and Rat. Recombinant immunogen, 207kDa. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	BTAf1 Rabbit pAb
Remarks/Alias	MOT1; TAF172; TAFII170; TAF(II)170; BTAf1
Species	Human
Gene ID (Human)	9044
Gene ID	9044
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA, ChIP
Cross Reactivity	Human, Rat
SWISS	O14981
Protein Weight	207kDa
Shipping	Ice bag