

## Anti-Btbd11 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Gqw0

Item Number: CM0022344



### Introduction

Rabbit polyclonal antibody targeting Btbd11, validated for WB and ELISA applications. Reacts with human, mouse, and rat. Btbd11 is an inhibitory interneuron-specific synaptic protein. Ideal for neuroscience research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Btbd11 Rabbit pAb                                                                                                 |
| Remarks/Alias    | Btbd11; 6330404E16Rik                                                                                             |
| Species          | Mouse                                                                                                             |
| Gene ID          | 74007 (Mouse)                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of mouse Btbd11 (NP_082985.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q6GQW0                                                                                                            |
| Protein Weight   | 122kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |