


Anti-C-Fos Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P01100

Item Number: CM0016528



Introduction

c-Fos rabbit polyclonal antibody validated for WB, IP, ELISA in human samples. Targets the AP-1 transcription factor subunit c-Fos (41 kDa), a key regulator of cell proliferation and differentiation.

[Learn More](#)

Parameter	Value
Product Name	c-Fos Rabbit pAb
Remarks/Alias	p55; AP-1; C-FOS; c-Fos
Species	Human
Gene ID (Human)	2353
Gene ID	2353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human c-Fos (NP_005243.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P01100
Protein Weight	41 kDa
Shipping	Ice bag