

Anti-C-Maf Polyclonal Antibody For Wb, Elisa - 075444

Item Number: CM0016767

Introduction

Rabbit polyclonal antibody against c-Maf (MAF) for WB and ELISA. Detects human, mouse, and rat c-Maf. Transcription factor involved in lens development, Th2 cell regulation, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	c-Maf Rabbit pAb
Remarks/Alias	CCA4; AYGRP; c-MAF; CTRCT21; c-Maf
Species	Human
GeneID	4094 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 314-403 of human c-Maf (NP_005351.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075444
Protein Weight	38kDa
Shipping	Ice bag