


Anti-C10Orf2/Twinkle Mtdna Helicase Polyclonal Antibody For Wb, Elisa - Q96Rr1

Item Number: CM0015766



Introduction

Rabbit polyclonal antibody targeting human C10orf2 (Twinkle mtDNA helicase), validated for Western blot and ELISA, cross-reacts with mouse. Suitable for mitochondrial DNA replication and repair research.

[Learn More](#)

Parameter	Value
Product Name	C10orf2 Rabbit pAb
Remarks/Alias	PEO; PEO1; SCA8; ATXN8; IOSCA; PEOA3; SANDO; TWINL; MTDPS7; PRLT55; C10orf2
Species	Human
Gene ID (Human)	56652
Gene ID	56652
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 585-684 of human C10orf2 (NP_068602.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96RR1
Protein Weight	77kDa
Shipping	Ice bag