


Anti-Ca3 Polyclonal Antibody For Wb, If/Icc, Elisa - P07451

Item Number: CM0015027



Introduction

High-quality anti-CA3 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat carbonic anhydrase III (CA3/CAIII), a muscle-specific cytoplasmic enzyme for reversible CO₂ hydration.

[Learn More](#)

Parameter	Value
Product Name	CA3 Rabbit pAb
Remarks/Alias	Car3, CAIII, CA3
Species	Human
GeneID (Human)	761
GeneID	761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CA3 (NP_005172.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07451
Protein Weight	30kDa
Shipping	Ice bag