

Anti-Cacna1S Polyclonal Antibody For If-P And Elisa - Q13698

Item Number: CM0019559



Introduction

Rabbit polyclonal antibody against human CACNA1S (Cav1.1, L-type calcium channel subunit alpha-1S) validated for IF-P and ELISA; cross-reacts with mouse and rat. Ideal for excitation-contraction coupling research.

[Learn More](#)

Parameter	Value
Product Name	CACNA1S Rabbit pAb
Remarks/Alias	MHS5; DHPRM; HOKPP; TTPP1; CMYP18; Cav1.1; HOKPP1; hypoPP; CCHL1A3; CACNL1A3; CACNA1S
Species	Human
GeneID (human)	779
GeneID	779
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1000 of human CACNA1S (NP_000060.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13698
Protein Weight	212kDa
Shipping	Ice bag