


Anti-Cadm3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N126

Item Number: CM0026663



Introduction

Rabbit polyclonal antibody against human CADM3 (NECL1, SynCAM3) validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience research of cell adhesion and synaptic organization.

[Learn More](#)

Parameter	Value
Product Name	CADM3 Rabbit pAb
Remarks/Alias	BigR; NECL1; TSLL1; CMT2FF; IGSF4B; Necl-1; synCAM3; CADM3
Species	Human
GeneID (Human)	57863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-330 of human CADM3 (NP_001120645.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N126
Protein Weight	43kDa
Shipping	Ice bag