

## Anti-Cald1 Rabbit Polyclonal Antibody For Wb And Elisa - Q05682

Item Number: CM0030963

### Introduction

Anti-CALD1 rabbit polyclonal antibody for WB and ELISA applications. Recognizes human caldesmon, an actin- and myosin-binding protein regulating smooth muscle contraction and non-muscle cell motility. Ideal for cardiovascular and oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CALD1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | CDM, HCAD, LCAD, h-CD, H-CAD, L-CAD, NAG22, CALD1                                                                |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 800                                                                                                              |
| Gene ID          | 800                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human CALD1 (NP_004333.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | Q05682                                                                                                           |
| Protein Weight   | 93kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |