


Anti-Camk2A Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Uqm7

Item Number: CM0006816



Introduction

Rabbit monoclonal antibody against human CAMK2A (CaMKII alpha, 54kDa) for Western blot, IF-paraffin, ELISA. Cross-reacts with mouse and rat. Ideal for studying synaptic plasticity, calcium signaling, and kinase activity.

[Learn More](#)

Parameter	Value
Product Name	CAMK2A Rabbit mAb
Remarks/Alias	CAMKA; MRD53; MRT63; CaMKIIalpha; CaMKIINalpha; CAMK2A
Species	Human
GeneID (Human)	815
Immunogen	Synthetic peptide within amino acids 201-300 of human CAMK2A (NP_741960.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQM7
Protein Weight	54kDa
Shipping	Ice bag