


Anti-Camki Rabbit Monoclonal Antibody For Wb And Elisa - Q14012

Item Number: CM0004890



Introduction

CaMKI Rabbit mAb validated for WB and ELISA with reactivity to human, mouse, and rat. Ideal for calcium/calmodulin-dependent signaling, neuronal development, and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	CaMKI Rabbit mAb
Remarks/Alias	CAMKI; CaMKI
Species	Human
GeneID (Human)	8536
GeneID	8536
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CaMKI (NP_003647.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14012
Protein Weight	41kDa
Shipping	Ice bag