

Anti-Camlg Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49069

Item Number: CM0015760



Introduction

Rabbit polyclonal antibody targeting human CAMLG (CAML/GET2), a key mediator of tail-anchored protein insertion and calcium signaling. Suitable for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	CAMLG Rabbit pAb
Remarks/Alias	CAML; GET2; CDG2Z; CAMLG
Species	Human
Gene ID (Human)	819
Gene ID	819
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human CAMLG (NP_001736.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49069
Protein Weight	33kDa
Shipping	Ice bag