


Anti-Caprin1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q14444

Item Number: CM0025549



Introduction

Rabbit polyclonal antibody against CAPRIN1, validated for Western blot, immunofluorescence/ICC, and ELISA. Reacts with human, mouse, and rat samples. Targets the RNA-binding protein involved in stress granule formation and mRNA regulation.

[Learn More](#)

Parameter	Value
Product Name	CAPRIN1 Rabbit pAb
Remarks/Alias	M11S1; GPIAP1; RNG105; GPIP137; GRIP137; p137GPI; CAPRIN1
Species	Human
Gene ID	4076
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CAPRIN1 (NP_005889.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14444
Protein Weight	78kDa
Shipping	Ice bag