


Anti-Carbonic Anhydrase 1 (Ca1) Monoclonal Antibody For Wb, Elisa - P00915

Item Number: CM0010086



Introduction

Rabbit monoclonal antibody against human Carbonic Anhydrase 1 (CA1), validated for WB and ELISA, cross-reactive with mouse and rat. Ideal for studying CO₂ metabolism and pH regulation.

[Learn More](#)

Parameter	Value
Product Name	Carbonic Anhydrase 1 (CA1) Rabbit mAb
Remarks/Alias	CAB; CA-I; Car1; HEL-S-11; Carbonic Anhydrase 1 (CA1)
Species	Human
GeneID (Human)	759
GeneID	759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-261 of human Carbonic Anhydrase 1 (CA1) (NP_001729.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P00915
Protein Weight	29kDa
Shipping	Ice bag